

¹²³I Iodine-123 Sodium Radiochemical Solution

Iodine-123 can be used for a wide range of clinical imaging applications in endocrinology, neurology, oncology and cardiology. It is also an ideal research tool because of its short half-life, clear SPECT images, and ease of handling.

BWXT Medical offers Sodium Iodide in solution and dry forms, with multiple shipment days available per week and next day delivery across North America.


Process Data

Nuclear Reaction

Primary



Secondary



Assay

High resolution gamma-ray spectrometry at 159 keV

Product Specification

Half-Life:

13.2 hours

Solution Form

Sodium iodide in 0.1N NaOH

Dry Form

Sodium iodide in 0.1N NaOH evaporated to dryness under nitrogen.

Specific Activity*

No carrier added

Activity Concentration*

At customer request

Radiochemical Purity

≥ 95% iodide

Radiopurity*

≥ 99.8% ¹²³I
≤ 0.2% total other nuclides

*At noon, Pacific Time, one day after shipment.

Please note that the radiochemicals listed are not tested for sterility or pyrogenicity. Verification of their suitability for use in humans is the sole responsibility of the purchaser.

Packaging

- V-vial with (1, 3, 5 mL)
- Glass vial (6, 10 mL)

Product Origin

Vancouver, Canada

Advantages



- Key supplier to North American market since 1984
- High quality product, preservative free
- Custom dispensing available for next day delivery
- Available in solution and dry forms



World-Class Delivery Service

BWXT Medical has a world-class distribution system and specialized order management system. Our extensive experience and carefully managed transportation logistics mean products reach your destination on time.

Other BWXT Products:

- Indium (In-111) Oxyquinoline Solution (US Only)
- Indium-111 Chloride Radiochemical Solution
- Iodine-123 Sodium Iodide Oral Solution (Canada Only)
- Germanium-68 Chloride Radiochemical Solution
- Strontium-82 Chloride Radiochemical Solution

Iodine-123 Decay Table Half-Life: 13.2 hours

HOURS	0	0.25	0.5	0.75	1	1.25	1.5	1.75	2	2.25	2.5	2.75
0	1.0000	0.9870	0.9741	0.9614	0.9488	0.9365	0.9243	0.9122	0.9003	0.8886	0.8770	0.8655
3	0.8542	0.8431	0.8321	0.8213	0.8105	0.8000	0.7895	0.7792	0.7691	0.7591	0.7492	0.7394
6	0.7297	0.7202	0.7108	0.7016	0.6924	0.6834	0.6745	0.6657	0.6570	0.6484	0.6400	0.6316
9	0.6234	0.6152	0.6072	0.5993	0.5915	0.5838	0.5762	0.5686	0.5612	0.5539	0.5467	0.5396
12	0.5325	0.5256	0.5187	0.5120	0.5053	0.4987	0.4922	0.4858	0.4794	0.4732	0.4670	0.4609
15	0.4549	0.4490	0.4431	0.4373	0.4316	0.4260	0.4204	0.4150	0.4096	0.4042	0.3989	0.3937
18	0.3886	0.3835	0.3785	0.3736	0.3687	0.3639	0.3592	0.3545	0.3499	0.3453	0.3408	0.3363
21	0.3320	0.3276	0.3234	0.3191	0.3150	0.3109	0.3068	0.3028	0.2989	0.2950	0.2911	0.2873
24	0.2836	0.2799	0.2762	0.2726	0.2691	0.2656	0.2621	0.2587	0.2553	0.2520	0.2487	0.2454
27	0.2422	0.2391	0.2360	0.2329	0.2299	0.2269	0.2239	0.2210	0.2181	0.2152	0.2124	0.2097
30	0.2069	0.2042	0.2016	0.1989	0.1964	0.1938	0.1913	0.1888	0.1863	0.1839	0.1815	0.1791
33	0.1768	0.1745	0.1722	0.1699	0.1677	0.1655	0.1634	0.1613	0.1592	0.1571	0.1550	0.1530
36	0.1510	0.1490	0.1471	0.1452	0.1433	0.1414	0.1396	0.1378	0.1360	0.1342	0.1324	0.1307
39	0.1290	0.1273	0.1257	0.1240	0.1224	0.1208	0.1192	0.1177	0.1161	0.1146	0.1131	0.1117
42	0.1102	0.1088	0.1073	0.1059	0.1046	0.1032	0.1019	0.1005	0.0992	0.0979	0.0966	0.0954
45	0.0941	0.0929	0.0917	0.0905	0.0893	0.0882	0.0870	0.0859	0.0848	0.0836	0.0826	0.0815
48	0.0804	0.0794	0.0783	0.0773	0.0736	0.0753	0.0743	0.0734	0.0724	0.0715	0.0705	0.0696
51	0.0687	0.0678	0.0669	0.0660	0.0652	0.0643	0.0635	0.0627	0.0618	0.0610	0.0602	0.0595
54	0.0587	0.0579	0.0572	0.0564	0.0557	0.0550	0.0542	0.0535	0.0528	0.0521	0.0515	0.0508
57	0.0501	0.0495	0.0488	0.0482	0.0479	0.0469	0.0463	0.0457	0.0451	0.0445	0.0440	0.0434

The half-life values stated are as of time of publication. BWXT values may vary.

People Strong

INNOVATION DRIVEN >

BWXT Medical Ltd. provides its customers, who conduct life-saving medical procedures for patients around the world, the benefit of decades of experience in the development, processing, packaging and delivery of medical isotopes, pharmaceuticals and medical devices. Headquartered in Kanata, Ontario, BWXT Medical employs approximately 300 highly skilled people in Kanata and Vancouver, British Columbia. BWXT Medical is part of the commercial operations segment of BWX Technologies, Inc. (NYSE:BWXT), headquartered in Lynchburg, Virginia, U.S.A. BWXT is a leading supplier of nuclear components and fuel to the U.S. government; provides technical and management services to support the U.S. government in the operation of complex facilities and environmental remediation activities; and supplies precision manufactured components, services and fuel for the commercial nuclear power industry. With approximately 6,600 employees, BWXT has 12 major operating sites in the U.S. and Canada. In addition, BWXT joint ventures provide management and operations at more than a dozen U.S. Department of Energy and NASA facilities. Follow us on Twitter @BWXTech and learn more at www.bwxt.com.

The information contained herein is provided for general information purposes only and is not intended nor to be construed as a warranty, an offer, or any representation of contractual or other legal responsibility.

The products and services described herein are provided by the subsidiaries of BWX Technologies, Inc.

© 2022 BWX Technologies, Inc. All rights reserved.



BWXT Medical
 447 March Road
 Ottawa, ON, Canada K2K 1X8

Regional Office
 4004 Wesbrook Mall
 Vancouver, BC, Canada V6T 2A3

Customer Service
isotopeorders@bwxt.com
 North America: +1.800.267.6211
 International: +1.613.963.2621
 General Emergency: +1.613.963.2300